


## YASARGIL VESSEL CLIP; STANDARD

permanent; 5mm; opening 7.2mm

200 gms; curved



**SKU:** 244-417-05

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE830

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.65

**SPECIFICATIONS** [curved](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE830

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.65

## SPECIFICATIONS

[curved](#), [permanent](#)