


YASARGIL VESSEL CLIP; STANDARD

permanent; 11mm; opening 8mm

200 gms; angled



SKU: 244-421-11

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
Medicon	58.55.41
SPECIFICATIONS	angled , permanent

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
Medicon	58.55.41
SPECIFICATIONS	angled , permanent