


YASARGIL VESSEL CLIP; STANDARD

permanent; 15mm; opening 9.2mm

200 gms; straight



SKU: 244-400-15

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE780

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE780

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)