


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 6.2mm

200 gms; curved



SKU: 244-423-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.42

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.42

SPECIFICATIONS [curved](#), [permanent](#)