


## YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 7.4mm

150 gms; straight; fenestrated Ø 3.5mm



**SKU:** 244-430-07

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE600
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.53.50
<b>SPECIFICATIONS</b>	<a href="#">permanent</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE600
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.53.50

## SPECIFICATIONS

[permanent](#), [straight](#)