


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 8.4mm

150 gms; straight; fenestrated Ø 3.5mm




SKU: 244-430-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE610
AUTHORS	YASARGIL
Medicon	58.53.51
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE610
AUTHORS	YASARGIL
Medicon	58.53.51

SPECIFICATIONS

[permanent](#), [straight](#)