


YASARGIL VESSEL CLIP; STANDARD

permanent; 12mm; opening 9.6mm

180 gms; straight; fenestrated Ø 3.5mm



SKU: 244-430-12

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

Aesculap FE620

AUTHORS [YASARGIL](#)

Medicon 58.53.52

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE620

AUTHORS [YASARGIL](#)

Medicon 58.53.52

SPECIFICATIONS

[permanent](#), [straight](#)