


YASARGIL VESSEL CLIP; STANDARD

permanent; 5mm; opening 5.8mm

150 gms; angled; fenestrated Ø 3.5mm



SKU: 244-431-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.53

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.53

SPECIFICATIONS [angled](#), [permanent](#)