


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 6.5mm

150 gms; angled; fenestrated Ø 3.5mm



SKU: 244-431-07

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.54

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.54

SPECIFICATIONS [angled](#), [permanent](#)