


YASARGIL VESSEL CLIP; STANDARD

permanent; 11mm; opening 7.2mm

180 gms: angled; fenestrated Ø 3.5mm



SKU: 244-431-11

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.55

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.55

SPECIFICATIONS [angled](#), [permanent](#)