


YASARGIL VESSEL CLIP; STANDARD

permanent; 17mm; opening 10.6mm

200 gms; straight




SKU: 244-400-17

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	permanent , straight