


## YASARGIL VESSEL CLIP; STANDARD

permanent; 17mm; opening 10.6mm

200 gms; straight



**SKU:** 244-400-17

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">permanent</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">permanent</a> , <a href="#">straight</a>