


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 7.9mm

150 gms; straight; fenestrated Ø 5mm



SKU: 244-436-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE640
AUTHORS	YASARGIL
Medicon	58.53.60
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE640
AUTHORS	YASARGIL
Medicon	58.53.60

SPECIFICATIONS

[permanent](#), [straight](#)