


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 9.1mm

150 gms; straight; fenestrated Ø 5mm



SKU: 244-436-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

Aesculap FE650

AUTHORS [YASARGIL](#)

Medicon 58.53.61

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE650

AUTHORS [YASARGIL](#)

Medicon 58.53.61

SPECIFICATIONS

[permanent](#), [straight](#)