


YASARGIL VESSEL CLIP; STANDARD

permanent; 12mm; opening 10.3mm

180 gms; straight; fenestrated Ø 5mm



SKU: 244-436-12

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE660

AUTHORS [YASARGIL](#)

Medicon 58.53.62

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE660

AUTHORS [YASARGIL](#)

Medicon 58.53.62

SPECIFICATIONS

[permanent](#), [straight](#)