


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 7.2mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-437-07

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.64

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.64

SPECIFICATIONS [angled](#), [permanent](#)