


## YASARGIL VESSEL CLIP; STANDARD

permanent; 11mm; opening 7.8mm

180 gms: angled; fenestrated Ø 5mm



**SKU:** 244-437-11

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.65

**SPECIFICATIONS** [angled](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.65

**SPECIFICATIONS** [angled](#), [permanent](#)