


YASARGIL VESSEL CLIP; STANDARD

permanent; 20mm; opening 11.4mm

200 gms; straight




SKU: 244-400-20

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	permanent , straight