


YASARGIL VESSEL CLIP; STANDARD

permanent; 11mm; opening 10mm

180 gms: angled; fenestrated Ø 5mm



SKU: 244-440-11

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.71

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.71

SPECIFICATIONS [angled](#), [permanent](#)