


YASARGIL VESSEL CLIP; STANDARD

permanent; 5mm; opening 5.8mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-441-05

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.74

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.74

SPECIFICATIONS [angled](#), [permanent](#)