


## YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 5.8mm

150 gms; angled; fenestrated Ø 5mm



**SKU:** 244-441-07

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN


**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.75

**SPECIFICATIONS** [angled](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.53.75

**SPECIFICATIONS** [angled](#), [permanent](#)