


YASARGIL VESSEL CLIP; STANDARD

permanent; 10mm; opening 5.8mm

150 gms; angled; fenestrated Ø 5mm



SKU: 244-443-10

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.79

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.53.79

SPECIFICATIONS [angled](#), [permanent](#)