


YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 6.2mm

110 gms; straight



SKU: 244-500-07

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.54.58

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.54.58

SPECIFICATIONS [straight](#), [temporary](#)