


YASARGIL VESSEL CLIP; STANDARD

temporary; 15mm; opening 9.2mm

90 gms; straight



SKU: 244-500-15

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE781

AUTHORS [YASARGIL](#)

Medicon 58.54.66

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE781

AUTHORS [YASARGIL](#)

Medicon 58.54.66

SPECIFICATIONS

[straight](#), [temporary](#)