

## YASARGIL VESSEL CLIP; STANDARD

temporary; 15mm; opening 9.2mm

90 gms; straight



**SKU:** 244-500-15

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE781
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.54.66
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">temporary</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE781
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.54.66

## SPECIFICATIONS

[straight](#), [temporary](#)