

## YASARGIL VESSEL CLIP; STANDARD

temporary; 17mm; opening 10.6mm

90 gms; straight



**SKU:** 244-500-17

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.54.67

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.54.67

**SPECIFICATIONS** [straight](#), [temporary](#)