

YASARGIL VESSEL CLIP; STANDARD

temporary; 25mm; opening 13.3mm

110 gms; straight



SKU: 244-500-25

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.54.72

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.54.72

SPECIFICATIONS [straight](#), [temporary](#)