


YASARGIL VESSEL CLIP; STANDARD

temporary; 9mm; opening 6.8mm

90 gms; curved



SKU: 244-501-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE753

AUTHORS [YASARGIL](#)

Medicon 58.55.60

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE753

AUTHORS [YASARGIL](#)

Medicon 58.55.60

SPECIFICATIONS

[curved](#), [temporary](#)