


YASARGIL VESSEL CLIP; STANDARD

temporary; 15mm; opening 8.7mm

110 gms; curved



SKU: 244-501-15

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.66

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.66

SPECIFICATIONS [curved](#), [temporary](#)