


YASARGIL VESSEL CLIP; STANDARD

temporary; 15mm; opening 7.4mm

110 gms; curved



SKU: 244-503-15

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.86

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.86

SPECIFICATIONS [curved](#), [temporary](#)