


## YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 5.4mm

110 gms; curved



**SKU:** 244-505-07

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.92

**SPECIFICATIONS** [curved](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.92

**SPECIFICATIONS** [curved](#), [temporary](#)