


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 8mm

110 gms; angled



**SKU:** 244-510-05T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Vessel Clips](#)

**Aesculap** FT270T

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FT270T

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

## SPECIFICATIONS

[angled](#), [temporary](#)