


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 8.7mm

90 gms; angled



SKU: 244-510-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)


Aesculap FT258T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.09

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT258T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.63.09

SPECIFICATIONS

[angled](#), [temporary](#)