


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 12mm; opening 10.2mm

90 gms; angled




SKU: 244-510-12T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.63.12
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.63.12

SPECIFICATIONS

[angled](#), [temporary](#)