


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 17mm; opening 11.9mm

90 gms; angled



**SKU:** 244-510-17T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.63.17
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">temporary</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>MATERIAL</b>	<a href="#">Titanium</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.63.17

## SPECIFICATIONS

[angled](#), [temporary](#)