


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 6.2mm

110 gms; angled



SKU: 244-511-05T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT103T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.85

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT103T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.61.85

SPECIFICATIONS

[angled](#), [temporary](#)