


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 6.2mm

110 gms; angled




SKU: 244-511-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT103T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.85
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT103T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.61.85

SPECIFICATIONS

[angled](#), [temporary](#)