


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 11mm; opening 8.4mm

110 gms; angled



SKU: 244-511-11T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.90
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.90

SPECIFICATIONS

[angled](#), [temporary](#)