


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 15mm; opening 8.7mm

200 gms; curved




SKU: 244-401-15T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.16
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.16

SPECIFICATIONS

[curved](#), [permanent](#)