


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 8mm; opening 5.4mm

130 gms; curved



**SKU:** 244-515-08T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.65.08

**SPECIFICATIONS** [curved](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.65.08

## SPECIFICATIONS

[curved](#), [temporary](#)