


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 7.2mm

90 gms; curved




SKU: 244-519-09T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.65.29 |
| SPECIFICATIONS | curved , temporary |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.65.29 |

SPECIFICATIONS

[curved](#), [temporary](#)