


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 20mm; opening 11mm

200 gms; curved




SKU: 244-401-20T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.20
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.61.20

SPECIFICATIONS

[curved](#), [permanent](#)