


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 9mm; opening 8.4mm

110 gms; straight; fenestrated Ø 3.5mm



**SKU:** 244-530-09T

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.64.09

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.64.09

## SPECIFICATIONS

[straight](#), [temporary](#)