


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 5.8mm

110 gms; angled; fenestrated Ø 3.5mm



SKU: 244-531-05T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.64.15
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.64.15

SPECIFICATIONS

[angled](#), [temporary](#)