


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 11mm; opening 7.2mm

110 gms; angled; fenestrated \varnothing 3.5mm



SKU: 244-531-11T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.21

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.21

SPECIFICATIONS

[angled](#), [temporary](#)