


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 4.9mm

110 gms; angled; fenestrated Ø 3.5mm



SKU: 244-533-07T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.27

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.27

SPECIFICATIONS

[angled](#), [temporary](#)