


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

permanent; 7mm; opening 6.2mm

150 gms; curved



SKU: 244-403-07T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT744T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.61.28

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT744T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.61.28

SPECIFICATIONS

[curved](#), [permanent](#)