


## YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 7mm; opening 7.9mm

110 gms; straight; fenestrated Ø 5mm



**SKU:** 244-536-07T

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)


**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.64.37

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**MATERIAL** [Titanium](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.64.37

## SPECIFICATIONS

[straight](#), [temporary](#)