


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 12mm; opening 10.3mm

90 gms; straight; fenestrated \varnothing 5mm



SKU: 244-536-12T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.42

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.64.42

SPECIFICATIONS

[straight](#), [temporary](#)