


YASARGIL VESSEL CLIP; STANDARD; TITANIUM

temporary; 5mm; opening 7.9mm

90 gms; angled; fenestrated Ø 5mm



SKU: 244-540-05T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|----------------------------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.64.65 |
| SPECIFICATIONS | angled , temporary |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.64.65 |

SPECIFICATIONS

[angled](#), [temporary](#)